ANKOM TDF Dietary Fiber Analyzer





- Automates AOAC 991.43, 985.29, 2001.03
- Optionally 2009.01, 2011.25, 2017.16 & 2022.01
- IDF/SDF/TDF Values
- Automates addition of chemicals, enzymes & rinses
- Auto controlled temperatures and agitation
- Reduces labor >50%
- Eliminates water baths and filtration flasks Frees up lab space
- 6 sample positions
- Accomplish 3 runs or 9 duplicates/day/shift (i.e. 18 samples within an 8 hr workday)

AOAC 991.43 – Conventional Steps

Prepare for Test		
1	Prepare Crucibles	
1a	Ash crucibles overnight	
1b	Clean with vacuum	
1c	Soak 1 hr @ room temp	
1d	Rinse with 3 solutions	
1e	Prepare DE filter bed	
1f	Dry at 130°C	
1g	Cool 1 Hour & weigh	
2	Weigh Samples	
3	Prepare/Heat Water Baths	
4	Place samples in beakers	

Digest Sample		
5	Heat to 95°C	
6	Add buffer	
7	Add Amylase	
8	Digest at 95°C for 35 min	
9	Scrape beaker walls and rinse with H ₂ O to move all residue to bottom	
10	Cool to 60°C	
11	Add Protease	
12	Digest at 60°C for 30 min	
13	Add HCI	
14	Add AMG	
15	Check/adjust pH	
16	Digest at 60°C for 30 min	

Isolate IDF Residue		
Transfer to filter		
Filter with vacuum		
Scratch filter bed as needed		
Rinse with 70°C H ₂ O		
Rinse with 78% EtOH		
Rinse with 95% EtOH		
Rinse with Acetone		
Dry IDF residue		
Weigh IDF residue		

Isolate SDF Residue		
26	Add 95% EtOH to filtrate	
27	Precipitate for 60 min	
28	Transfer to filter	
29	Filter with vacuum	
30	Scratch filter bed as needed	
31	Rinse with 78% EtOH	
32	Rinse with 95% EtOH	
33	Rinse with Acetone	
34	Dry SDF residue	
35	Weigh SDF residue	

Calculate Results		
36	Correct for Protein	
37	Correct for Ash	
38	Calculate results using weights and corrections	

AOAC 991.43 – Ankom Automated

Prepare for Test		
1	Weigh Filter Bags	
2	Weigh DE and Samples	
3	Prepare/Heat Water Baths	
4	Place samples in beakers	

Digest Sample		
5	Heat to 95°C	
6	Add buffer	
7	Add Amylase	
8	Digest at 95°C for 35 min	
10	Cool to 60°C	
11	Add Protease	
12	Digest at 60°C for 30 min	
13	Add HCI	
14	Add AMG	
15	Check/adjust pH	
16	Digest at 60°C for 30 min	

Isolate IDF Residue		
17	Transfer to filter	
18	Filter with vacuum	
20	Rinse with 70°C H ₂ O	
21	Rinse with 78% EtOH	
22	Rinse with 95% EtOH	
23	Rinse with Acetone	
24	Dry IDF residue	
25	Weigh IDF residue	

Isolate SDF Residue			
26	Add 95% EtOH to filtrate		
27	Precipitate for 60 min		
28	Transfer to filter		
29	Filter with vacuum		
31	Rinse with 78% EtOH		
32	Rinse with 95% EtOH		
33	Rinse with Acetone		
34	Dry SDF residue		
35	Weigh SDF residue		

Calculate Results		
36	Correct for Protein	
37	Correct for Ash	
38	Calculate results using weights and corrections	

Blank = Steps eliminated due to automation and Filter Bag Technology Green = Steps automated by instrument, no technician intervention White = Steps still required by technician

AOAC 991.43 – Ankom Automated

Prepare for Test		Dig	est Sample
1	Weigh Filter Bags		
2	Weigh DE and Samples		
3	Install Filter Bags		
4	Place sample & DE in bags		
		15	Check/adjust pH

Isolate IDF Residue		
23	Rinse with Acetone	
24	Dry IDF residue	
25	Weigh IDF residue	

Isolate SDF Residue

33	Rinse with Acetone
34	Dry SDF residue
35	Weigh SDF residue

Calculate Results	
36	Correct for Protein
37	Correct for Ash
38	Calculate results using weights and corrections

Went from 38 manual steps to 14

AOAC 2017.19 – Conventional Steps

Prepare for Test

1	Prepare Crucibles
1a	Ash crucibles overnight
1b	Clean with vacuum
1c	Soak 1 hr @ room temp
1d	Rinse with 3 solutions
1e	Prepare DE filter bed
1f	Dry at 130°C
1g	Cool 1 Hour & weigh
2	Weigh Samples
3	Prepare water baths
4	Place samples in beakers

Digest Sample		
5	Add buffer	
6	Add Amylase/AMG	
7	Heat to 37°C	
8	Digest for 4 hours	
9	Scrape beaker walls and	
	rinse with H ₂ O to move	
	all residue to bottom	
10	Add Trizma Base	
11	Check/adjust pH	
12	Heat to 90°C	
13	Incubate for 20 min	
14	Cool to 60°C	
15	Add Protease	
16	Digest for 30 min	
17	Add Acetic Acid	
18	Check/ adjust pH	
19	Add internal standard	

Isolate IDF Residue	
20	Transfer to filter
21	Filter with vacuum
22	Scratch filter bed as needed
23	Rinse with H ₂ O
24	Rinse with 78% EtOH
25	Rinse with 95% EtOH
26	Rinse with Acetone
27	Dry IDF residue
28	Weigh IDF residue

Isolate SDF Residue		
29	Add 95% EtOH to filtrate	
30	Precipitate for 60 min	
31	Transfer to filter	
32	Filter with vacuum	
33	Scratch filter bed as needed	
34	Rinse with 78% EtOH	
35	Rinse with 95% EtOH	
36	Rinse with Acetone	
37	Dry SDF residue	
38	Weigh SDF residue	

Calculate Results		
39	Correct for Protein	
40	Correct for Ash	
41	Determine WASDF	
42	Calculate results using weights and corrections	

AOAC 2017.19 – Ankom Automated

Prepare for Test	
1	Weigh filter bags
2	Weigh Samples
3	Prepare water baths
4	Place samples in beakers

Blank = Steps eliminated due to automation		
and Filter Bag Technology		
Green = Steps automated by instrument, no		
technician intervention		
White = Steps still required by technician		

Digest Sample	
5	Add buffer
6	Add Amylase/AMG
7	Heat to 37°C
8	Digest for 4 hours
10	Add Trizma Base
11	Check/adjust pH
12	Heat to 90°C
13	Incubate for 20 min
14	Cool to 60°C
15	Add Protease
16	Digest for 30 min
17	Add Acetic Acid
18	Check/ adjust pH
19	Add internal standard

Isolate IDF Residue		
20	Transfer to filter	
21	Filter with vacuum	
23	Rinse with H ₂ O	
24	Rinse with 78% EtOH	
25	Rinse with 95% EtOH	
26	Rinse with Acetone	
27	Dry IDF residue	
28	Weigh IDF residue	

Isolate SDF Residue	
29	Add 95% EtOH to filtrate
30	Precipitate for 60 min
31	Transfer to filter
32	Filter with vacuum
34	Rinse with 78% EtOH
35	Rinse with 95% EtOH
36	Rinse with Acetone
37	Dry SDF residue
38	Weigh SDF residue

Calcu	culate Results			
39	Correct for Protein			
40	Correct for Ash			
41	Determine WASDF			
42	Calculate results using weights and corrections			

AOAC 2017.19 – Ankom Automated

Prepare for Test		Digest Sample	
1	Weigh filter bags		
2	Weigh Samples	11	Check/adjust pH
3	Prepare water baths		
4	Place samples in beakers		
_		18	Check/ adjust pH

19

Add internal standard

Isolate IDF Residue		
20	Transfer to filter	
21	Filter with vacuum	
26	Rinse with Acetone	
27	Dry IDF residue	
28	Weigh IDF residue	

Isolate SDF Residue		
36	Rinse with Acetone	
37	37 Dry SDF residue	
38	Weigh SDF residue	

Calcu	alculate Results			
39	Correct for Protein			
40	Correct for Ash			
41	Determine WASDF			
42	Calculate results using weights and corrections			

Went from 42 manual steps to 19

ANKOM TDF Dietary Fiber Analyzer



Filter Bag Technology

- Dual chamber Filter Bag Technology
 - Allows the separation of Digestion and Filtration
 - Eliminates sample loss during transfer steps
- Increased filtering surface area reduces the time needed to filter samples and eliminates the use of vacuum flasks and crucibles
- Enables automation
- Reduces error
- Reduces labor and cost per assay

Insoluble Dietary Fiber (IDF) Filter Bag Dietary Fiber Flow-Thru Bag Soluble Dietary Fiber (SDF) Filter Bag



Comparison of Data – IDF



Comparison of Data – SDF



Comparison of Data – TDF

