

DVE™ Viscometer

Our most affordable digital viscometer

What's Included

- Instrument
 - 6 spindles (RV/HA/HB) or 4 spindles (LV)
 - Carrying Case
 - Lab Stand (Model G)
 - Spindle Guard Leg*
- *Not applicable to HA / HB torque models

Optional Accessories

- Small Sample Adapter
- UL Adapter
- Thermosel
- DIN Adapter
- Viscosity Standards
- RV/HA/HB-1 Spindle
- Quick Action Lab Stand
- Temperature Bath
- Quick Connect/Extension Links



Features

User-friendly Interface

- Keypad control
- Sharp viewing screen for close up or distance viewing

No Calculations Required

Direct reading of viscosity

Displayed Info:

- Viscosity (cP, P, mPa•s or Pa•s)
- % Torque
- Speed/Spindle

Easy-to-Use

Range:

Push for determining full scale range (FSR) viscosity

Bubble Level:

Conveniently located for easy adjustment

Accuracy: $\pm 1.0\%$ of range

Repeatability: $\pm 0.2\%$

18 Speeds

For complete range capability

VISCOSITY RANGE cP(mPa•s)			SPEEDS (200 Available)	
Model	Min	Max	RPM	Number of Increments
DVELV	1†	2M	.3-100	18
DVERV	100††	13M	.3-100	18
DVEHA	200††	26M	.3-100	18
DVEHB	800††	104M	.3-100	18

† 1 cP achieved with UL Adapter accessory.

15 cP on LV with standard spindles.

†† Minimum viscosity is achieved with optional RV/HA/HB-1 spindle.

M = 1 million cP = Centipoise

mPa•s = Millipascal•seconds

