

-40°C Freezer



Technical Parameters

Model	BDF-40V362	
External Size(W*D*H)	690*744*1885mm	
Internal Size(W*D*H)	520*567*1266mm	
Capacity	362L	
Type	Vertical	
Temp. Range	-15°C~-40°C	
Temp. Accuracy	0.1°C	
Control System	Microprocessor Control; LED Display	
Alarm	High or low temperature alarm, door ajar alarm, power failure alarm, sensor failure alarm, battery switch alarm	
Refrigeration Type	Direct cooling	
Refrigerant	R290	
Compressor	Well-known brand	
Construction	Structure	Uniboday design
	Internal Material External Material	Spraying steel plate
Shelf	4pcs	
Consumption	360W	
Power Supply	AC220V±10%, 50Hz; 110V±10%, 60Hz	
Package Size	810*860*2100mm	
Gross Weight	165kg	

-40°C Freezer



58L

368L

Technical Parameters

Model	BDF-40V58		BDF-40V368
External Size(W*D*H)	540*620*1115mm		720*830*1870mm
Internal Size(W*D*H)	361*405*406mm		Upper:580*626*625mm Down:466*636*625mm
Capacity	58L		Upper:189; Down:179L
Type	Vertical		
Temp. Range	-15°C~-40°C		Upper:2°C~8°C Down:-20°C~-40°C
Temp. Accuracy	0.1°C		
Control System	Microprocessor Control LED Display		
Alarm	BDF-40V58&BDF-40V398: High or low temperature alarm, sensor failure alarm; BDF-40V368: High or low temperature alarm, sensor failure alarm, power failure alarm, door open alarm.		
Refrigeration Type	Direct Refrigeration		
Refrigerant	R404a and R23		R600a and R290
Compressor	Unicompressor, imported famous brand compressor		
Material	Structure	Uniboday design	
	Internal	304 stainless steel	Sprayed steel plate
	External	Sprayed steel plate	Sprayed steel plate
Shelf	1 pc shelf		Upper:3 pcs shelves Down:2 pcs shelves and 6 drawers
Consumption	340W		402W
Power Supply	AC220V±10%, 50/60HZ		
Package Size	610*690*1275mm		775*920*1990mm
Gross Weight	98kg		166kg